

Name Index Number. /

233/3
CHEMISTRY
Paper 3
PRACTICAL
Oct./ Nov. 2012
2¼ hours

Candidate's Signature:

Date:



THE KENYA NATIONAL EXAMINATIONS COUNCIL
Kenya Certificate of Secondary Education
CHEMISTRY
Paper 3
PRACTICAL
2¼ hours

233/3 - Chemistry Paper - 3	
Wednesday	8.00 am - 10.15 am
31/10/2012	(1st Session)

Instructions to Candidates

- (a) Write your name and index number in the spaces provided above.
- (b) Sign and write the date of examination in the spaces provided above.
- (c) Answer **ALL** the questions in the spaces provided.
- (d) You are not allowed to start working with the apparatus for the first 15 minutes of the 2¼ hours allowed for this paper. This time is to enable you to read the question paper and make sure you have all the chemicals and apparatus that you may need.
- (e) Mathematical tables and silent electronic calculators may be used.
- (f) All working **MUST** be clearly shown where necessary.
- (g) This paper consists of 8 printed pages.
- (h) Candidates should check the question paper to ascertain that all the pages are printed as indicated and that no questions are missing.

For Examiner's use only

Question	Maximum Score	Candidate's Score
1	20	
2	11	
3	09	
Total Score	40	

© 2012 THE KENYA NATIONAL EXAMINATIONS COUNCIL
 Kenya Certificate of Secondary Education, 2012
 CHEMISTRY

Paper 3

1. You are provided with:

- solution A containing an oxidising agent A;
- solution B, 0.05 M aqueous sodium thiosulphate;
- solution C containing a reducing agent C;
- aqueous potassium iodide;
- solution D, starch solution.

You are required to determine the:

- concentration of solution A;
- rate of reaction between the oxidising agent A and the reducing agent C.

Procedure 1

1. Using a pipette and **pipette filler**, place 25.0 cm³ of solution A into a 250 ml conical flask.
2. Measure 10 cm³ of aqueous potassium iodide and add it to solution A in the conical flask. Shake the mixture. Add 10 cm³ of 2 M sulphuric (VI) acid to the mixture and shake.
3. * Fill a burette with solution B and use it to titrate the mixture in the conical flask until it just turns **orange-yellow**. Add 2 cm³ of solution D to the mixture in the conical flask. Shake thoroughly. Continue titrating until the mixture **just turns colourless**. Record your results in **table 1** below.
4. Repeat the procedure and complete table 1. **Retain the remainder of solution A and solution D** for use in procedure II.

Table 1	I	II	III
Final burette reading	16.30	16.30	16.30
Initial burette reading	0.00	0.00	0.00
Volume of solution B used (cm ³)	16.30	16.30	16.30

(4 marks)

(a) Calculate the:

- (i) average volume of solution B used; (1 mark)
-

- (ii) number of moles of sodium thiosulphate. (1 mark)
-

3

(b) Given that one mole of A reacts with six moles of sodium thiosulphate, calculate the:

(i) number of moles of A that were used; (1 mark)

.....

(ii) concentration of solution A in moles per litre. (2 marks)

.....

.....


Procedure II

1. Label six test - tubes as 1, 2, 3, 4, 5 and 6 and place them in a test - tube rack.
2. Using a clean burette, measure the volumes of distilled water shown in **table 2** into the labelled test - tubes.
3. Using a burette, measure the volumes of solution **A** shown in **table 2** into each of the test - tubes.
4. Clean the burette and rinse it with about 5 cm³ of solution **C**.
5. Using the burette, measure 5 cm³ of solution **C** and place it into a 100 ml beaker.
6. Using a 10 ml measuring cylinder, measure 5 cm³ of solution **D** and add it to the beaker containing solution **C**. Shake the mixture.
7. Pour the contents of test - tube number 1 to the mixture in the beaker and immediately start a stop watch. Swirl the contents of the beaker. Record the time taken for a **blue** colour to appear in **table 2**.
8. Repeat steps 5 to 7 using the contents of test - tube numbers 2, 3, 4, 5 and 6.
9. Complete **table 2** by computing $\text{Rate} = \frac{1}{\text{time}} (\text{S}^{-1})$.

Table 2

Test - tube number	1	2	3	4	5	6
Volume of distilled water (cm ³)	0	2	3	5	6	7
Volume of solution A (cm ³)	10	8	7	5	4	3
Time (seconds)						
Rate = $\frac{1}{\text{time}} (\text{S}^{-1})$						

(a) Plot a graph of rate (y-axis) against volume of solution A. (3 marks)



(b) What time would be taken for the blue colour to appear if the experiment was repeated using 4 cm³ of distilled water and 6 cm³ of solution A? (2 marks)

.....

.....

.....

2. You are provided with solid E. Carry out the experiments below. Write your observations and inferences in the spaces provided.

(a) Place all of solid E in a boiling tube. Add about 20 cm³ of distilled water and shake until all the solid dissolves, label the solution as solution E. Use solution E for experiments (i) and (ii).

(i) To 2 cm³ of solution E, in a test - tube in each of experiments I, II, III and IV, add:

I. two drops of aqueous sodium sulphate;

Observations

Inferences

(1 mark)

(1 mark)

II. five drops of aqueous sodium chloride;

Observations

Inferences

(1 mark)

(1 mark)

III. two drops of barium nitrate;

Observations

Inferences

(1 mark)

(1 mark)

IV. two drops of lead (II) nitrate;

Observations

Inferences

(1 mark)

(1 mark)

7

- (ii) To 2 cm³ of solution E, in a test - tube, add 5 drops of aqueous sodium hydroxide. Add the piece of aluminium foil provided to the mixture and shake. Warm the mixture and test any gas produced with both blue and red litmus papers.

Observations	Inferences
(2 marks)	(1 mark)

3. You are provided with solid F. Carry out the following tests. Write your observations and inferences in the spaces provided.

- (a) Place all of solid F in a boiling tube. Add about 20 cm³ of distilled water and shake until all the solid dissolves. Label the solution as solution F.

Add about half of the solid sodium hydrogen carbonate provided to 2 cm³ of solution F.

Observations	Inferences
(1 mark)	(1 mark)

- (b) (i) Add about 10 cm³ of dilute hydrochloric acid to the rest of solution F in the boiling tube. Filter the mixture. Wash the residue with about 2 cm³ of distilled water. Dry the residue between filter papers. Place about one third of the dry residue on a **metallic** spatula and burn it in a Bunsen burner flame.

Observations	Inferences
(1 mark)	(1 mark)

8

- (ii) Place all the remaining residue into a boiling tube. Add about 10 cm³ of distilled water and shake thoroughly. **Retain the mixture for the tests in (C).**

Observations	Inferences
(½ mark)	(½ mark)

- (c) Divide the mixture into two portions:

- (i) to the first portion, add the rest of the solid sodium hydrogen carbonate.

Observations	Inferences
(1 mark)	(1 mark)

- (ii) to the second portion, add two drops of bromine water.

Observations	Inferences
(1 mark)	(1 mark)

THIS IS THE LAST PRINTED PAGE.